



Features

- Keyed temperature lock out
- Push button digital operation
- Large dual digital displays
- Displays actual temperature and set point at all times
- Separate digitally controlled heating chambers
- Rapid warm time (2-6 hours depending on load)
- Uniform heating
- Fully insulated
- Audible and Visual Over-Temp Alarm
- U.L./C.U.L. approved
- FDA Registered
- OSHPD Pre-Approved



Typical
Manufacturer's
Image



Axia SWC 24 Blanket Warmer

Specifications

Overall Size

- 30"W x 26.5"D x 24.5"T
- 1 Adjustable Shelf

Chamber I.D.

- 26"W x 23"D x 15.25"T
- 5.27 cubic feet of storage
- Temperature range: 90-160 F.
- Accuracy +/- 1 degree F. of set point

Electrical

- 120V, 7.5 amp., 50/60HZ
- .30 kw/hr, 1024 btu/hr
- UL/CUL listed
- FDA registered

Seattle Technology: Surgical Division
3915 152nd St. NE Marysville, WA 98271
info@st surg.com
1-800-827-3747



Note: The technical data given in this publication is for general information and are subject to change without notice. Actual configuration on the unit may vary. Contact our sales representatives for a complete list of details.