



The Dräger Caleo Incubator provides an ideal microenvironment system for premature babies and sick infants up to a body weight of 5 kg (11 lbs.) or a body length of 55 cm (22 inches), providing a controlled environment of warmth, humidity, and elevated O<sub>2</sub> concentration in the patient area. The total body weight when treating twins is limited to 5 kg (11 lbs.). Caleo is intended for use in clinical environments where premature babies or infants are treated who require a controlled climatic environment.



# Dräger Caleo

## Options for nursing and therapy:

- Convective heat therapy through control of incubator air temperature or infant skin temperature
- Humidification of the incubator air
- O2 therapy through controlled elevation of the O2 concentration in the patient environment
- Nursery and intensive care via hand ports or two large access doors
- Pivoting bed for raising and lowering the infant's head (Trendelenburg and anti-Trendelenburg position)

## With monitoring for:

- Air temperature
- Skin temperature
- Relative humidity
- O2 concentration
- Weight

## Dimensions

Incubator (width x depth)	1167 mm x 687 mm (45.9" x 27.0")
Height overall with height adjustable column	1220 mm to 1520 mm (48.0" x 59.8")
Height overall with fixed column	1270 mm / 1370 mm / 1470 mm (50" / 53.9" / 57.9")
Height of mattress surface (variable height)	800 mm (31.5") to 1100 mm (43.3")
Height of mattress surface (fixed height)	875 mm (34.4") / 950 mm (37.4") / 1050 mm (41.3"), depending on fixed height column installed
Bed (width x depth)	645 mm x 500 mm (25.4" x 19.7")

## Weight

Overall weight max.	230 kg (507 lbs)
basic configuration	130 kg (286.6 lbs)
Overall load capacity	66 kg (142.5 lbs)

Seattle Technology: Surgical Division  
3915 152nd St. NE Marysville, WA 98271  
info@st surg.com  
1-800-827-3747



Note: The technical data given in this publication is for general information and is subject to change without notice. Actual configuration on the unit may vary. Contact our sales representatives for a complete list of details.