



GE OEC 8800 FlexiView C-Arm

Learning from the strength of thousands of systems installed worldwide and years of intraoperative imaging experience, the GE OEC 8800 FlexiView has created a flexible and economical mobile C-arm with the high quality imaging technology of our premium platform and the simplicity and maneuverability of a compact point-and-shoot system. Now you can experience superb image quality for your conventional intraoperative imaging needs.

- 1k x 1k x 16 bit image processing
- Touchscreen controls
- AutoTrak™ ABS control
- Removable digital image storage (static and dynamic)
- Instant laser quality film printer:
 - Integrated DICOM connectivity
 - Infra-red remote control

GE OEC 8800 FlexiView C-Arm

Maximum Focal Spot	1.5 mm
Motorized Tilting/Pivot Rotation	370
X-ray Generator Power Rating	2.2 kV



3915 152nd St. NE Marysville, WA 98271
1-800-827-3747
www.stsurg.com
info@stsurg.com