



Seattle Technology  
3915 152nd St NE  
Marysville, WA 98271

1.800.827.3747

stsurg.com

## Product Overview

The Axia SCW 24 single door blanket warmer is the perfect answer for warming blankets. Designed with optimal patient comfort in mind, this blanket warmer will fit well in most healthcare environments.

## Features

- Keyed temperature lock out
- Push button digital operation
- Large dual digital displays
- Displays actual temperature and set point at all times
- Separate digitally controlled heating chambers
- Rapid warm time (2-6 hours depending on load)
- Uniform heating
- Fully insulated
- Audible and Visual Over-Temp Alarm
- U.L./C.U.L. approved
- FDA Registered
- OSHPD Pre-Approved



## Specifications

### Overall Size

30"W x 26.5"D x 24.5"T  
1 Adjustable Shelf

### Chamber I.D.

26"W x 23"D x 15.25"T  
5.27 cubic feet of storage  
Temperature range: 90-160 F.  
Accuracy +/- 1 degree F. of set point

### Electrical

1120V, 7.5 amp., 50/60HZ  
.30 kw/hr, 1024 btu/hr  
UL/CUL listed  
FDA registered

Note: The technical data given in this publication is for general information and is subject to change without notice. Actual configuration on the unit may vary. Contact our sales representatives for a complete list of details.