

Product Overview

SIREMOBIL® Compact is equipped with a high-duty X-ray generator to meet even the most demanding clinical applications. Challenging projections, like lateral lumbar, which require penetration of dense anatomy, are easily mastered. When cardiology applications such as EP studies and pacemakers placement demand exceptionally long fluoro times, you're ready. SIREMOBIL® Compact's uniquely designed X-ray tube provides fluoro times virtually without limit – over 50 minutes – with high image quality and lowest possible dose. Long fluoro times are even more important as mobile C-arms are used increasingly as back-up systems for other modalities.



Siemens Siremobil Compact C-Arm

Features

- Versatile
- Compact size
- High image quality
- Single key operation
- Highly mobile
- Dual Mode 9/6 with Digistore 100 and laser targeting

Specifications on Next Page

Note: The technical data given in this publication is for general information and is subject to change without notice. Actual configuration on the unit may vary. Contact our sales representatives for a complete list of details.



Seattle Technology
3915 152nd St NE
Marysville, WA 98271
1.800.827.3747
stsurg.com

Specifications

Versatility

The capabilities and performance of SIREMOBIL® Compact have established it as a reliable, versatile & economical mobile C-arm designed to meet the full requirements of today's surgery. Whether used in traditional applications like orthopedics, or in the growing areas of vascular intervention, gastrointestinal surgery or minimally invasive procedures, SIREMOBIL® Compact provides advanced features, thoughtful design and Siemens' quality.

Compact Size

Transportability and maneuverability are optimized with its small footprint and the extended, Counterbalanced C-arm design. The system features a streamlined user interface for fast, efficient and easy operation.

Image Quality

The SIREMOBIL® brand name is known for its superior imaging components, and SIREMOBIL® Compact carries on the legacy. While image quality is maximized, all operating modes are optimized toward dose savings.

Subtraction and Roadmap

The vascular memory offers subtraction at the touch of a button for accurate visualization during interventional procedures, as well as roadmap for fast, easy positioning during catheterization and dilatation.

Networking

Keeping pace with the industry's development in networking, SIREMOBIL® Compact can be equipped with on and off-line communication capabilities. High Duty, Lower Dose.

Note: The technical data given in this publication is for general information and is subject to change without notice. Actual configuration on the unit may vary. Contact our sales representatives for a complete list of details.