

## Product Overview

Amico OR Pendants incorporate multiple equipment rails located the front and back of console to allow integration of multiple accessories. The OR Pendant is capable of withstanding up to 400 lbs with adjustable equipment shelf yielding maximum weight of 150 lb each. Integrated with the OR Pendant, our intuitive control handle allows for easy maneuvering of the equipment rack.



**Amico OR Pendant Equipment Boom**

## Features

- Adjustable shelf and drawer configurations available
- Shelves hold up to 150 lbs
- Customizable console with integrated rails, medgas and electrical outlets
- 100% customizable: including size as well as location of electrical and medgas services
- We recommend blank provisions to allow the addition of future service
- Mount accessories in any one of the integrated rails
- All accessories are modular and easily moved on the custom Amico rail
- Provides easy access to all internal components for hassle free maintenance
- Prevents dust and bacteria from collecting on the top surface of the console
- Easy to clean round cover
- No need to pipe in Instrument Air or Nitrogen – unlike pneumatic brakes
- Intuitive control handle for easy ergonomic use
- Energy efficient

Note: The technical data given in this publication is for general information and is subject to change without notice. Actual configuration on the unit may vary. Contact our sales representatives for a complete list of details.