

Product Overview

The AMSCO 3080 RL surgical table is one of the most durable, flexible and reliable surgical table available. It features adaptable patient position with simple level commands to increase efficiency in the operating room. This all-round surgery table combines easy patient positioning and good c-arm access to accommodate the most routine laparoscopic procedure to the most demanding orthopedic case to neurosurgery.

The Quantum 3080 RL also has a proven track record for outstanding reliability. Even in the event of a complete power failure, the user can still position the patient or return the table to the neutral position.



AMSCO Quantum 3080 RL Surgical Table

Features

- Self-leveling brakes
- One-press Return to Level
- X-ray cassette deck
- Pressure reduction mattress
- Manual override foot pump
- Optional foot control
- Optional battery back up
- Will accommodate the most routine laparoscopic procedure to the most demanding orthopedic case to neurosurgery
- Complimented by a radiolucent tabletop
- Removable head section to enhance available imaging window
- Absolute Override Protection
Even in the event of a complete power failure and electrical system malfunction on the table, the user can still position the patient or return the table to the neutral position. Every table is equipped with a manual foot pump and auxiliary switches that provide override control of all critical patient positions.

Specifications on Next Page



Seattle Technology
3915 152nd St NE
Marysville, WA 98271

1.800.827.3747

stsurg.com

Specifications

Raised Height:
44 inches

Lowered Height:
27 inches

Length:
75 15/16 inches

Width:
20 inches

Width with Rail:
21 5/8 inches

Weight Limit:
500 lbs

Trendelenberg:
25°

Reverse:
25°

Trendelenberg Reflex:
20°

Flex:
30°

Lateral Tilt:
±18°

Back:

Down
25°

Up
55°

Leg:

Down
105°

Up
80°

Note: The technical data given in this publication is for general information and is subject to change without notice. Actual configuration on the unit may vary. Contact our sales representatives for a complete list of details.