

Product Overview

The AMSCO 3080 SP Surgical Table is one of the most durable, flexible, and reliable surgical tables available. It features simple patient positioning with simple level commands to increase efficiency in the operating room. The 3080 SP also has good C-Arm access and convenience of use. Plus, the 3080 SP has a proven track record of reliability. Use of the hand control makes positioning the patient very intuitive. Using the Patient Orientation button makes positioning the table simple when reversing the patient on the table. When the procedure is completed a single depression of the Return to Level button brings the table back to the neutral (level) position.

Features

- Self-leveling brakes
- One-press Return to Level
- X-ray cassette deck
- Pressure reduction mattress
- Manual override foot pump
- Optional foot control
- Optional battery back up
- Will accommodate the most routine laparoscopic procedure to the most demanding orthopedic case to neurosurgery
- Complimented by a radiolucent tabletop
- Removable head section to enhance available imaging window
- Absolute Override Protection
Even in the event of a complete power failure and electrical system malfunction on the table, the user can still position the patient or return the table to the neutral position. Every table is equipped with a manual foot pump and auxiliary switches that provide override control of all critical patient positions.



AMSCO Quantum 3080 SP Surgical Table

Specifications on Next Page



Seattle Technology
3915 152nd St NE
Marysville, WA 98271

1.800.827.3747

stsurg.com

Specifications

Raised Height:
44 inches

Lowered Height:
27 inches

Length:
75 15/16 inches

Width:
20 inches

Width with Rail:
21 5/8 inches

Weight Limit:
500 lbs

Trendelenberg:
25°

Reverse:
25°

Trendelenberg Reflex:
20°

Flex:
30°

Lateral Tilt:
±18°

Back:

Down
25°

Up
55°

Leg:

Down
105°

Up
80°

Note: The technical data given in this publication is for general information and is subject to change without notice. Actual configuration on the unit may vary. Contact our sales representatives for a complete list of details.