



Seattle Technology
 3915 152nd St NE
 Marysville, WA 98271
 1.800.827.3747
 stsurg.com

Product Overview

The Chick 703 Orthopedic and Surgical Operating Table provides more for less...

- It's the only low price table with a full-length table top.
- The top provides body support for other surgical procedures and allows C-arm imaging of upper body extremities.
- Low acquisition cost and diversified surgical applications make the 703 extremely productive and profitable.
- It's the only table in its class with a patient transfer board as standard.
- The top rotates 180° for upper extremity traction.
- Its exclusive, retractable leg spars and light weight reduce storage space and physical effort.

The Chick 703. Your best investment when you want a top quality, full featured, versatile table at an economical price.

Standard Components:

- Radiolucent perineal post 2 inches.
- Cast saddle, adult.
- Safety strap hook and pile, white.
- Perineal post pad for 2 inch post.
- Disposable pad set (1 only).
- Cushion set (6 piece).
- Radiolucent head section.
- Shoulder board, radiolucent w/o cushion.
- Radiolucent lumbar section.
- Saddle adapter.
- Well leg support assembly.
- Right Bell foot piece w/traction unit.
- Radiolucent saddle assembly.
- Left bell foot piece w/traction unit.
- Stabilizer kit.
- Patient transfer board.
- Radiolucent counter traction support/horizontal perineal post w/pads.



Chick 703 Orthopedic & Surgical Operating Table

Standard Components Continued:

- Tibia extension bracket.
- Safety strap, side rail mount, black.
- 90° pin and wire holder
- Side rail extension, pair.

Specifications:

Table Height:
 Raised..... 46 inches
 Lowered..... 36 inches

Table Top:

Length with radiolucent head section, shoulder section and head section:..... 51.5 inches
 Width:..... 20 inches

Table Width:..... 28 inches
 Table Length:..... 62-123 inches

Traction Bars:..... 11-34 inches

Table Top Indexing: 180° w/Lock

Leg Spars:
 Adduction..... 0-60°
 Abduction..... 0-74°

Brakes:..... Manual

Weight Limit:..... 350 lbs

Note: The technical data given in this publication is for general information and is subject to change without notice. Actual configuration on the unit may vary. Contact our sales representatives for a complete list of details.