



Seattle Technology
 3915 152nd St NE
 Marysville, WA 98271
 1.800.827.3747
 stsurg.com

Product Overview

The Maquet Alphastar has unprecedented flexibility when it comes to patient positioning. The head and leg sections are totally interchangeable, not to mention the wide variety of additional table sections that are available to optimize patient positioning.

Among the many attachments are a carbon fiber extension that provides un-impaired visibility when using a C-arm. The leg section can also be replaced with traction bars making the table a good choice for trauma centers.

The strong column and wide footprint allow the patient to be positioned in normal and reversed orientation with out compromising stability. Combine this with the flexibility and connectivity of the patient platform and this table provides the ability to position any patient optimally for surgery without compromising patient comfort.

This table also offers a totally unique head support that has been designed specifically to meet the needs of Maxillo-Facial, Ophthalmology and ENT procedures.

The Alphastar provides the ability to transport the table (complete with patient when lowered and leveled) using the Auto Drive feature. The drive system requires no pushing or pulling by the care giver to move the table and patient, thereby reducing the risk of backstrain.



Maquet Alphastar Plus Surgery Table

Specifications:

Table Height:
 Raised..... 42 inches
 Lowered (Ultra)..... 27 inches
 Lowered (Lo Version)... 23 1/2 inches

Trendelenberg:..... 30°
 Reverse Trendelenberg:..... 30°

Head Section:
 Raise..... 0-15°
 Lower..... 0-90°

Lateral Tilt:..... ±18°

Back Section:
 Down..... 25°
 Up..... 68°

Seat Section:
 Flex..... 60°
 Reflex..... 70°

Leg Section (Removable):
 Down..... 105°
 Up..... 80°

Return to Level..... '0' Button

Auto Drive..... Hand Control

Brakes (Electronic) Hand Control

Weight Limit:..... 500 lbs

Length:..... 82 1/2 inches
 Width:..... 20 1/2 inches
 Width w/Rail:..... 23 1/2 inches

Note: The technical data given in this publication is for general information and is subject to change without notice. Actual configuration on the unit may vary. Contact our sales representatives for a complete list of details.