



Seattle Technology
 3915 152nd St NE
 Marysville, WA 98271
 1.800.827.3747
 stsurg.com

Product Overview

The Skytron 6500 Elite General Purpose Surgical Table further extends the range of applications of the 6000 series. It increases the patient weight limit to 500 lbs as well as increasing the degrees of articulation offered. It also incorporates a powered Kidney Bridge as standard along with fully powered Brakes and a "Return-to-Level" function, all of which are controlled from the Hand Control.

The rotating top of 210° further enhances C-arm and microscope access along with more leg room for a seated surgeon and in many cases better surgical access. A stable footprint and rigid surface provides assurance that the patient position will be maintained as desired.

The table can also be upgraded to add an x-ray top providing unobstructed passage of x-ray film cassettes up to 14 inches x 17 inches.

The 6500 Model also features a detachable head section that can be inserted into the end of the foot section to provide support for the legs when positioning patients for lithotomy procedures.

This model is also equipped with a manual brake release that is independent of the main controls and requires no power to operate.



SKYTRON 6500 Elite General Purpose Surgical Table

Specifications:

Table Height:		Seat Section:	
Raised.....	47 inches	Flex.....	0-70°
Lowered.....	26 inches	Reflex.....	0-90°
Trendelenberg:.....	30°	Kidney Bridge:	5 inch Elevation
Reverse Trendelenberg:....	30°	Leg Section:	
Head Section:		Down.....	105°
Raise.....	0-15°	Up.....	20°
Lower.....	0-90°	Table Top Rotation:.....	210°
Lateral Tilt:.....	±30°	Weight Limit:.....	500 lbs
Back Section:		Length:.....	76 inches
Down.....	40°	Width:.....	20 inches
Up.....	90°	Width w/Rail:	21 3/4 inches
		Battery:.....	Optional

Note: The technical data given in this publication is for general information and is subject to change without notice. Actual configuration on the unit may vary. Contact our sales representatives for a complete list of details.