



Seattle Technology
 3915 152nd St NE
 Marysville, WA 98271
 1.800.827.3747
 stsurg.com

Table Overview

The Skytron Hercules 6700B General Purpose Surgery Table provides unparalleled imaging access for today's most demanding cases while lifting 1200 lbs and articulating an unprecedented 1000 lbs. The removable Back and Leg Sections facilitate the use of a Carbon Fiber Back and Leg Section to further extend the imaging window.

Due to the limited mobility of the patients being catered for by this table, as well as the challenges faced by caregivers during patient transfer this table features an "ultra-low" Low Position of 23 inches.

The rotating top of 210° further enhances C-arm and microscope access along with more leg room for a seated surgeon and in many cases better surgical access. A stable footprint and rigid surface provides assurance that the patient position will be maintained as desired.

The table can also be upgraded to add an x-ray top providing unobstructed passage of x-ray film cassettes up to 14 inches x 17 inches.

The Hercules 6700B table has been fully tested to 4X its rated capacity to provide caregivers with total confidence when positioning extremely heavy patients.

This model also offers the convenience of an optional wireless IR Remote Control.

Finally, for additional stability the 6700 Series features a 4-point "auto-leveling" braking system with the brakes being located at the outer corners of the base.



SKYTRON 6700B Hercules Surgical Table

Specifications:

Table Height:
 Raised..... 44 inches
 Lowered..... 23 inches

Trendelenberg:..... 30°
 Reverse Trendelenberg:..... 30°

Head Section:
 Raise..... 0-15°
 Lower..... 0-90°

Lateral Tilt:..... ±30°

Back Section:
 Down..... 40°
 Up..... 90°

Seat Section:
 Flex..... 0-40°
 Reflex..... 0-90°

Kidney Bridge:.... 6 inch Elevation

Leg Section:
 Down..... 100°
 Up..... Horizontal

Table Top Rotation:..... 210°

Weight Limit (Lift):..... 1200 lbs
 Weight Limit (Articulate): 1000 lbs

Length:..... 76 inches
 Width:..... 20 inches
 Width w/Rail:..... 22 inches

Battery:..... Optional

Return to Level:..... Yes

Note: The technical data given in this publication is for general information and is subject to change without notice. Actual configuration on the unit may vary. Contact our sales representatives for a complete list of details.