


Product Overview

BiCap III Electrosurgical Unit (ESU) incorporates ConMed's Dynamic Response technology, which reacts to different tissue impedances in less than 1 millisecond ensuring a consistent delivery of power and clinical effect while minimizing surrounding tissue damage.



Features

- ESUm mono and bipolar
- 200W power
- Dynamic response technology reacts to different tissue impedances ensuring a consistent delivery of power and clinical effect while minimizing surrounding tissue damage.
- Front panel volume control
- Option to select automatic programming which restores the BiCap III to last settings used
- Pulse Cut Mode provides precise energy bursts for critical dissections

Specifications on Next Page

Note: The technical data given in this publication is for general information and is subject to change without notice. Actual configuration on the unit may vary. Contact our sales representatives for a complete list of details.



Seattle Technology
3915 152nd St NE
Marysville, WA 98271

1.800.827.3747

stsurg.com

Specifications

BiCap® III Monopolar Footswitch, Low Profile

Cable Length 15.00ft (4.57m)

BiCap® III ESU 120.00V

Modes

Pure cut

Coag

Blend cut

Bipolar

Pulse cut

Optional

BiCap® III System Including Footpedals

Note: The technical data given in this publication is for general information and is subject to change without notice. Actual configuration on the unit may vary. Contact our sales representatives for a complete list of details.