

Product Overview

Philips IntelliVue MP5SC with IntelliVue Guardian Early Warning Scoring (EWS) helps enhance patient care by providing sought-after measurements and Clinical Decision Support right at the point of care.



Philips IntelliVue MP5SC Spot Check Monitor

Features

- **Essential measurements to support informed decisions**

Philips is committed to providing excellent standard measurements such as NBP, Microstream CO₂ (including integrated Pulmonary Index), choice of spot check temperature technology with Suretemp, Tympanic Temperature or Temporal Temperature and a choice of pulse oximetry with Philips FAST SpO₂, Nellcor™ OxiMax™ SpO₂ or Masimo® rainbow® SET® SpO₂ (including certain Masimo rainbow parameters).

- **Flexible operational modes for use on the general ward**

The monitor can be used in spot-check or continuous monitoring mode. In continuous mode, the IntelliVue MP5SC can be used for long-term bedside monitoring. Its optional EWS Escalation package allows the MP5SC to convert to a full ECG monitor via a patient-worn device when enhanced monitoring is required, such as during a code event.

Features Continued on Next Page



Seattle Technology
3915 152nd St NE
Marysville, WA 98271
1.800.827.3747
stsurg.com

Features Continued

- **Large 8.4" touchscreen shows key data at a glance**

This monitor has a large and bright 8.4" touchscreen with flexible parameter configurations to enhance patient care.

- **Flexible MX40 connection with short-range radio**

Connects to the IntelliVue MX40 via Short Range Radio for multi-patient vital signs and occasional viewing of continuous monitoring.

- **IntelliVue user interface simplifies monitoring**

The IntelliVue user interface is optimized for ease of use on the general floor. It supports efficient workflow, helping caregivers in a variety of settings to simply measure, evaluate, and document their patient's vital signs and EWS scores.

- **Standard supplies – same as rest of IntelliVue monitors**

The MP5SC uses the same accessories as the rest of the IntelliVue family, standardizing training, ordering, and use throughout your facility.

- **IntelliVue Guardian Early Warning Scoring to aid rapid response**

The IntelliVue MP5SC is fully configurable with IntelliVue Guardian Early Warning Scoring (EWS) settings, parameters, alarms, and screen layouts. It provides caregivers with direct access to an automated scoring system based on their own hospital's policy and guidelines. Whatever the level of patient acuity and frequency of vital sign collection, IntelliVue Guardian EWS can help caregivers identify patient deterioration early with fast and easy Early Warning Scoring calculations.

- **Clinical Decision Support – Standard IntelliVue Guardian EWS**

Standard IntelliVue Guardian EWS Clinical Decision Support with optional IntelliVue Guardian EWS Escalation Package for patient-worn, device-based wireless ECG monitoring.

Note: The technical data given in this publication is for general information and is subject to change without notice. Actual configuration on the unit may vary. Contact our sales representatives for a complete list of details.