

Product Overview

Axia's UroMax 4 C-arm surgical tables offer an excellent solution for real-time urological X-ray procedures. The cantilevered, carbon fiber, radiolucent top areas allow full fluoroscopic visualization and unobstructed C-arm positioning. The UroMax 4 has highly versatile movements, options, and accessories. The carbon fiber leg extension accessory can quickly turn your Urology table into a procedure table for pain management and other C-arm guided procedures.



Axia UroMax 4 Surgical Table

Features

- Carbon fiber with integral facial cutout section
- Leg extension easy to remove and very light weight
- Diving board style for easy C-arm access
- Large patient capacity (500-lb pound load)
- A 2" patient high-quality comfort pad
- Light weight – easy to physically move the table
- Our 24" wide top eliminates the needs for arm boards
- Easy patient transfer to table
- Low height for patient boarding
- ***Stirrups are an optional accessory***

Specifications

Dimensions:

Table Length:
55 in (139 cm)

Table Width:
24 in (60.9 cm)

Table Height:
30 – 40 in (76.2 – 101.6 cm)

Weight Capacity:
500 lbs (226.7 kg)

Specifications Continued on Next Page



Seattle Technology
3915 152nd St NE
Marysville, WA 98271
1.800.827.3747
stsurg.com

Specifications Continued

Table Weight:
740 lbs (335.7 kg)

H-Shaped Steel Base:
51 x 26 in (129.5 – 66 cm)

Mattress:
2 in thick (5 cm)

Radiolucent Area:
27 x 24 in (68.5 x 60.9 cm)

Table Top:
55 x 24 x 1.62 in (139 x 60 x 4 cm)

Tilts:

Lateral Tilt:
 $\pm 17^\circ$

Trendelenburg:
 $\pm 13^\circ$

Longitudinal Travel:
10 in (± 5 in)

Power:

AC Voltage:
120V~

Frequency:
60 Hz

AC Current:
6.5A

Power Continued:

Electric Shock Protection:
Class I

Applied Part Rating:
Type B

Mode of Operation:
10% Duty Cycle; 1 min on, 9 min off

Environmental:

Storage Humidity:
0- 92%

Operating Humidity:
0- 92%

Storage Temperature:
32°F – 104°F (0°C – 40°C)

Operating Temperature:
Up to 104°F (40°C)

Note: The technical data given in this publication is for general information and is subject to change without notice. Actual configuration on the unit may vary. Contact our sales representatives for a complete list of details.