

Product Overview

The Tuttnauer 2340M Autoclave will satisfy all of your sterilization needs without compromising on quality, safety or reliability. The Tuttnauer 2340M manual autoclave is designed to compliment any healthcare facility. The bright new easy to use panel was designed with the operator in mind.



Tuttnauer 2340M Autoclave

Features

- Automatic shut off at the end of both the sterilization and dry cycles.
- A long life electro polished chamber and door.
- A double safety locking device prevents the door from opening while chamber is pressurized.
- Drain valve is located on the front, allowing for quick and easy draining of water reservoir.
- Dual safety thermostat to protect against overheating.

International certifications.

Specifications on Next Page

Note: The technical data given in this publication is for general information and is subject to change without notice. Actual configuration on the unit may vary. Contact our sales representatives for a complete list of details.



Seattle Technology
3915 152nd St NE
Marysville, WA 98271

1.800.827.3747

stsurg.com

Specifications

Chamber Dimensions:

x L – inches (mm) 9" x 18 (230 x 470)

Chamber Volume:

5 Gal – 19 Liters

Overall Dimensions DxWxH – inches (mm):

21.5" x 20" x 14.4" – (545 x 510 x 365)

Standard Cassettes Capacity:

2 full and 2 half

Tray Dimensions – inches (mm):

16.3" x 6.7" x 0.8" – (415 x 170 x 20)

No. of Trays:

3

Standard Unwrapped Cycle Time:

Cold: 27 Min. Hot: 13 Min.

Voltage:

(V) – Frequency (Hz) 120 – 50/60

Power:

(W) – Current (A) 1400W – 12A

Autoclave Weight:

– Lbs (kg) 85 (39)

Compliance:

PED 97/23 – MDD 93/42 – ASME – UL/EN
61010-1 – FDA – ISO 13485:2003 – ISO
9001:2000

Note: The technical data given in this publication is for general information and is subject to change without notice. Actual configuration on the unit may vary. Contact our sales representatives for a complete list of details.