



Seattle Technology  
3915 152nd St NE  
Marysville, WA 98271

1.800.827.3747

stsurg.com

## Product Overview

The Getinge 46-Series Washer has been specially developed to match the needs of larger surgical instrument sets in central sterile supply departments. It cleans and disinfects all kinds of items on five levels. Each wash level, consisting of a removable shelf and spray arm, can be adapted to the individual need. Tubular instruments can be cleaned in the MIS wash cart with capacity for as many as 80 connections.



## Features

- Unique drying unit that is effective, fast, and energy efficient
- Ergonomic handle
- Glass door makes it easy to overview the process
- Easy to clean, there are no external screws that can gather dirt
- AN-wash cart with capacity for 4–6 sets
- Slim design featuring a width of only 25.5 inches
- Large Chambers
- A washer-disinfector with suitable accessories to meet your demands

***Specifications on the Next Page***

Note: The technical data given in this publication is for general information and is subject to change without notice. Actual configuration on the unit may vary. Contact our sales representatives for a complete list of details.



Seattle Technology  
3915 152nd St NE  
Marysville, WA 98271

1.800.827.3747

stsurg.com

## Specifications

	<b>46-4</b>	<b>46-5</b>
● <b>External Dimensions</b>		
Without Dryer (HxWxD)	63 x 25.5 x 27 inches	66 x 25.5 x 27 inches
With Dryer (HxWxD)	73 x 25.5 x 27 inches	76 x 25.5 x 27 inches
Chamber Dimensions (HxWxD)	23 x 22 x 24 inches	26 x 22 x 24 inches
Chamber Volume	9.88 cubic feet	10.77 cubic feet
Effective Chamber Volume	7.1 cubic feet	7.94 cubic feet
DIN 1/1 Trays per Load	8	10
Wash Levels	4	5
Weight	Without Dryer 331 lbs With Dryer 397 lbs	Without Dryer 342 lbs With Dryer 406 lbs
● <b>Water</b>		
Water Supply	CW and WW standard, TW option	
Consumption per phase	5-6 gallons (depending on the wash-cart)	
Detergent Pumps	2 standard (+2 optional)	
Circulation Pump Capacity	1.1 kW	

Note: The technical data given in this publication is for general information and is subject to change without notice. Actual configuration on the unit may vary. Contact our sales representatives for a complete list of details.