

Product Overview

These 18 DIN trays capacity machines are the largest Steelco washer disinfector models designed to meet the increased reprocessing needs of CSSD. These machines assure consistently tested and efficient washing-disinfecting performances. These models are available in standard or Fast Cycle configurations, saving cycle time (total cycle time reduced up to 30 minutes) and reaching higher levels of energy and water savings. The ergonomic design of the door level height allows the user a convenient loading/unloading job. A powerful vertically installed washing pump ensures a high water flow rate combined with effective spray pressure and also grants a complete drain of the washing circuit.



Features

- Washing arms Low friction bearings ease chamber and cart washing arms rotation for an improved efficiency of the washing and drying phases.
- Washing cart automatic identification.
- Triple water filtering system captures residue preventing re-circulation.
- Active water coupling system The connection system of the hydraulic circuit, between the washer and the washing cart, allows a tight seal junction for a more efficient cleaning and uniform disinfection temperature coverage.

Specifications on the Next Page

Note: The technical data given in this publication is for general information and is subject to change without notice. Actual configuration on the unit may vary. Contact our sales representatives for a complete list of details.



Seattle Technology
3915 152nd St NE
Marysville, WA 98271
1.800.827.3747
stsurg.com

Specifications

- **Models:**

DS 900 – Single door washer-disinfector

DS 1000 – Double door washer-disinfector

DS 900 1S / DS 900 2S / DS 900 3S
Single door, fast cycle washer disinfectors

DS 1000 1S / DS 1000 2S / DS 1000 3S
Double door, fast cycle washer disinfectors

- **Washer Dimensions**

W D H: mm 1100 960 1940

inches 43.30 37.80 76.38

Same dimensions for the Fast Cycle versions.

- **Load Capacity**

18 DIN Trays

- **Chamber Volume**

500 lt / ~132.09 Gal US

- **Wash Cart Volume**

350 lt / ~92.46 Gal US

- **Narrow Frame Configuration**

All models can also be configured into a narrow frame to reduce their width to 900mm only.

Please check with your Regional Sales Manager for availability.

- **Washer Dimensions for Standard Cycle**

(WxDxH): mm 900 x 960 x 1940 / inches 35.43 x 37.80 x 76.38

- **Washer Dimensions for Fast Cycle**

(WxDxH): mm 900 x 960 x 2620 / inches 35.43 x 37.80 x 104.72

Note: The technical data given in this publication is for general information and is subject to change without notice. Actual configuration on the unit may vary. Contact our sales representatives for a complete list of details.