

Product Overview

The 3003 Impulse puts automatic control of the imaging surface right at your fingertips.

Featuring joystick control, the table offers 65" of unobstructed imaging. The two-speed, multi-directional top slide takes the work out of manual maneuvers and automatically manipulates the imaging surface.



Skytron 3003 Impulse Surgical Table

Features

- 20" longitudinal and 8" lateral top slide
- 500 lb lift and articulation weight capacity
- Two-speed, multidirectional joystick control
- 28"- 42" adjustable table height
- 83" carbon fiber top
- 65" unobstructed imaging window
- 20° lateral tilt
- 20° Trendelenburg and reverse Trendelenburg

Specifications on Next Page

Note: The technical data given in this publication is for general information and is subject to change without notice. Actual configuration on the unit may vary. Contact our sales representatives for a complete list of details.



Seattle Technology
3915 152nd St NE
Marysville, WA 98271

1.800.827.3747

stsurg.com

Specifications

Dimensions

Length: 83" (2,100 mm)

Width: 19 5/8" (500 mm)

Width (With Side Rails): 22" (560 mm)

Adjustments Via Hand Control

Height: 28"-42" 710-1065 mm

Lateral tilt (left/right): 20°

Trendelenburg: 20°

Reverse Trendelenburg: 20°

Joystick Control

For operation, the top button on joystick control must be depressed

Orientation Button: Reverse Normal

Center Button: Moves table top to center position

Speed: High Low

Max Permissible Load

500 lb (227kg) Lift

500 lb (227kg) Articulation

Note: The technical data given in this publication is for general information and is subject to change without notice. Actual configuration on the unit may vary. Contact our sales representatives for a complete list of details.