

Product Overview

The Aaron Smoke Shark II is designed to be easy to use and extremely effective at smoke and particle removal at the operative site. While designed for use with all electrosurgical generators, it is also excellent for use in controlling laser smoke.

The Bovie Aaron Smoke Shark II Smoke Evacuator has a modern, light-weight, compact design allowing for simplified use, storage and stack-ability. The SF35 Extended-Life Filter offers up to 35 hours of performance reducing procedure cost by ~38% from the previous model when used at the lowest setting. The covered 3-port filter design accommodates a variety of tubing sizes [1/4"(6.4mm), 3/8"(9.5mm), and 7/8"(22mm)] providing versatility in the procedure and accessory setup. The ultra-quiet design offers the same excellent performance with less noise.



Bovie Aaron Smoke Shark II Smoke Evacuator

A variety of accessories are available to customize the Smoke Shark to your specific application. These include laser-resistant tubing, tubing reducers for attachment to smoke pencils or vaginal speculum, many other tubing configurations (sterile or nonsterile), foot control and a remote switch activator.

Features

- A modern, light-weight, compact design allowing for simplified use, storage and stackability
- Never again order the Smoke Shark and still need to order a filter upon receipt, the SE02 is packaged with both the SF35 Extended-Life Filter and FSSEP Pneumatic Footswitch.
- The SF35 Extended-Life Filter offers up to 35 hours of performance reducing procedure cost by 38%(+/-) from the previous model, when used at the lowest setting.

Features Continued on Next Page



Seattle Technology
3915 152nd St NE
Marysville, WA 98271
1.800.827.3747
stsurg.com

Features

- Covered 3-port filter design accommodates a variety of tubing sizes [1/4"(6.4mm), 3/8"(9.5mm), and 7/8"(22mm)] providing versatility in procedure and accessory setup
- The ultra-quiet design offers the same excellent performance with less noise
- The SE02 Smoke Shark Smoke Evacuator is -3 and RoHS compliant
- **Complete PACKAGE**
 - Unit
 - Power Cord
 - User's Guide
 - Pneumatic Footswitch (FSSEP)
 - Long Life 35-hour Filter (SF35)

Specifications

Line Voltage

SE02: 100/120 VAC

SE02-220: 220/240 VAC

Line Frequency

50-60 Hz

Dimensions and Weight

Width: 27.9 cm (11")

Height: 15.2 cm (6")

Depth: 39.4 cm (15.5")

Weight: 4.3 kg (10 lbs.)

Note: The technical data given in this publication is for general information and is subject to change without notice. Actual configuration on the unit may vary. Contact our sales representatives for a complete list of details.