



Seattle Technology
3915 152nd St NE
Marysville, WA 98271

1.800.827.3747

stsurg.com

Product Overview

The Schure Spine Frame provides a stable and convenient platform for spine surgeries such as laminectomies, decompressions, and disc surgery. You can adjust it intraoperatively to the desired lordosis to get optimal surgical site exposure. Tapered adjustable pads provide optimal patient safety and comfort. The Schure Spine Frame has improved features over the current Wilson™ Frame, has a patient weight capacity of 500 lbs. and is C-arm compatible. It can be used with any O.R. table, and fits the Jackson Spine Table and Allen Flex Frame as well. Comes complete with transport storage dolly and crank handle, disposable headrest and 2 arm cradle pads.



Schure Spine Frame

Features

- Adjusts intra-operatively to the desired lordosis to get optimal surgical site exposure
- 500 lb patient weight capacity
- Can be used with any O.R. table and fits the Jackson Spine Table and Allen Flex Frame as well
- Frame comes complete with Schure™ Spine Frame Dolly, Crank Handle, Disposable Head Rest, and
- 2-Arm Cradles
- Radiolucent and allows for unrestricted C-arm access
- Deluxe foam pads included

Specifications on Next Page

Note: The technical data given in this publication is for general information and is subject to change without notice. Actual configuration on the unit may vary. Contact our sales representatives for a complete list of details.



Seattle Technology
3915 152nd St NE
Marysville, WA 98271

1.800.827.3747

stsurg.com

Specifications

Dimensions (W x D x H):

18" x 7" x 30"

Weight:

18 lbs



Note: The technical data given in this publication is for general information and is subject to change without notice. Actual configuration on the unit may vary. Contact our sales representatives for a complete list of details.