

Product Overview

The Steris 4085 Surgical Table is a mobile, electro-hydraulically operated surgical table specifically designed to provide patient positioning flexibility required for modern surgical care facilities. The narrow column, radiolucent table, and slide provide a large image footprint. It has outstanding access for the C-arm and surgical team without time-consuming patient reorientation. The imaging area permits clear, high-quality images that are crucial in today's MIS, Cardiothoracic, orthopedic, and neurosurgical procedures.



AMSCO 4085 Surgical Table

Features

- Accommodates virtually all patients with generous weight and height ranges
- Provides outstanding access for C-arm, permitting high quality images crucial to cardiothoracic, orthopedic and neurosurgical procedures
- Contours to virtually all patient postures with four section tabletop
- Radiolucent 18" sliding top for total body imaging (43" upper and 43" lower body C-arm access)
- Provides unrestricted perineal access with removable leg section
- Enhances lateral procedures with a powered, radiolucent kidney elevator
- Easy to use, easy to read hand control with digital data window

Specifications on Next Page

Note: The technical data given in this publication is for general information and is subject to change without notice. Actual configuration on the unit may vary. Contact our sales representatives for a complete list of details.



Seattle Technology
3915 152nd St NE
Marysville, WA 98271

1.800.827.3747

stsurg.com

Specifications

General

Table Length 81" (2057 mm)

Table Width 20" (508 mm) (28" with extenders)

Table Height Range 26" to 45" (660 to 1143 mm)

Table Slide Range 9" (227 mm) to head, 9" to foot
(18" total)

Table Shipping Weight 560 lbs. (254 kg)

Patient Weight Capacity

1,100 lbs. (500 kg) patient support, including
raise/lower (centered on the column)

1,000 lbs (454 kg) full table articulation (centered on
the column)

600 lbs (272 kg) full table articulation, including slide

Tilting Features

Trendelenburg / Reverse 30° / 30°

Lateral Tilt (left/right) 20° / 20°

Head Section +90° to -90°

Back Section +80° to -40°

Leg Section 0/-105° (removable)

Flex/Reflex 140° / 100°

Manual Override Yes

Perneal Cut-out Yes

Kidney Elevator 4" powered radiolucent

Note: The technical data given in this publication is for general information and is subject to change without notice. Actual configuration on the unit may vary. Contact our sales representatives for a complete list of details.